

## Dermojet hr

PRODUCT REFERENCE: 474



The Dermojet type HR replaces the standard Dermojet for medical use with the advantage that it can be sterilised at 134° C for 18 minutes.

The Dermojet type HR is equipped with:

- nozzle body made of stainless steel, a central piston, a tension spring and a tension lever
- an injection head
- a Pyrex or Radel® reservoir (capacity 4 ml)
- a stainless steel tip

Weight: approximately 300 g

Dosage: Set to a single 0.1 ml dose.

The Dermojet HR has the same general characteristics and uses as the standard Dermojet except that the stainless steel tip and Pyrex reservoir replace the polysulfone used in the standard Dermojet which cannot withstand autoclave sterilisation.

## **Applications**

- Local anaesthesia
- Hair transplant
- Corticoid injection
- The Dermojet HR works without cartridge, bottle or ampoule
- The Dermojet HR does not come with any medicinal substances; it is the practitioner who chooses the solution to be injected when filling the reservoir.